


Hero Flex Order Form

Specification

Wrist Type



☐ Falcon Wrist



☐ USMC Wrist

Chirality


☐ Left

☐ Right

Affected Arm

For guidance on any of these measurements, please refer to our Clinician Guidelines, or contact our configuration team at configuration@openbionics.com.

A - Olecranon to Residuum End (cm)



Sound Arm

E - Olecranon to Terminal Device (cm)

This is the distance between the elbow bone and the base of the terminal device (end of Hero-Flex).



There must be a minimum difference in length between measurements A and E of 2.8mm, or select 'Make as short as possible'.

OR

Make as short as possible

Additional Options

Pull-through hole?

☐ Yes ☐ No

Please note, depending on the patient's residual limb, we may need to add additional length to accommodate a pull-through. We will notify you during the design process if this is the case.

Add Kidney Condyles?

☐ Yes ☐ No


Kidney condyles attach to the inside of the frame between the frame and the socket, to aid suspension. *This is not available for the short version of the Hero Flex.

MyoBand?

☐ Yes ☐ No

Add a MyoBand to the order - this allows you to place MyoPods outside of the socket for forearm suspended Flex.

Notes - Please include any important details regarding the arm configuration here:

A large, light gray rounded rectangular box with a thin black border, intended for the user to provide additional details regarding the arm configuration.

Hero Flex Cast/Scan Information

Cast Type

Please let us know if you'll be sending a scan file, or a physical cast. Please note, we cannot work from a check socket or direct scan of a patient's limb.

☐ Digital Scan

☐ Physical Cast

Cast/Scan Notes - Please include any important details regarding your cast/scan here:

Please contact your sales rep for the address to ship your cast to.

PLEASE NOTE: Casts will be destroyed after scanning. If you require the cast to be returned, please indicate this below. You will need to arrange a collection once your cast has been scanned.

Cast handling preference:

☐ Destroy cast after scanning

☐ Return cast (I will arrange and cover collection costs)