


# Hero Flex Order Form

## Specification

### Wrist Type



☐ Falcon Wrist



☐ USMC Wrist

### Chirality


☐ Left

☐ Right

### Affected Arm

For guidance on any of these measurements, please refer to our Clinician Guidelines, or contact our configuration team at [configuration@openbionics.com](mailto:configuration@openbionics.com).

A - Olecranon to Residuum End (cm)



### Sound Arm

E - Olecranon to Terminal Device (cm)

This is the distance between the elbow bone and the base of the terminal device (end of Hero-Flex).



**There must be a minimum difference in length between measurements A and E of 2.8mm, or select 'Make as short as possible'.**

OR

Make as short as possible

## Additional Options

### Pull-through hole?

☐ Yes ☐ No

Please note, depending on the patient's residual limb, we may need to add additional length to accommodate a pull-through. We will notify you during the design process if this is the case.

### Add Kidney Condyles?

☐ Yes ☐ No


Kidney condyles attach to the inside of the frame between the frame and the socket, to aid suspension. \*This is not available for the short version of the Hero Flex.

### MyoBand?

☐ Yes ☐ No

Add a MyoBand to the order - this allows you to place MyoPods outside of the socket for forearm suspended Flex.

Notes - Please include any important details regarding the arm configuration here:

A large, light gray rounded rectangular box with a thin black border, intended for the user to provide notes regarding the arm configuration.

# Hero Flex Cast/Scan Information

## Cast Type

Please let us know if you'll be sending a scan file, or a physical cast. Please note, we cannot work from a check socket or direct scan of a patient's limb.

☐ Digital Scan

☐ Physical Cast

**Cast/Scan Notes** - Please include any important details regarding your cast/scan here:

**Please contact your sales rep for the address to ship your cast to.**

**PLEASE NOTE:** Casts will be destroyed after scanning. If you require the cast to be returned, please indicate this below. You will need to arrange a collection once your cast has been scanned.

Cast handling preference:

☐ Destroy cast after scanning

☐ Return cast (I will arrange and cover collection costs)

# Hero RGD/PRO Order Form

## Hand Selection:

☐ Hero PRO      ☐ Hero RGD

## Chirality

☐ Left      ☐

Right

## Size

☐ Medium      ☐ Large

☐ Add a Socket Component Kit?

**Please only select this option if you are not also ordering a Hero Flex and would like to make your own socket.**

Includes : 1 x Falcon wrist, 2 x Falcon lamination ring, 4 x MyoRays (for inserting MyoPods into socket) and lamination cap.