



MAINTENANCE AND CARE

- Brace should be worn over a cotton shirt
- Fold the lumbar sacral brace and attach the handle near the power units when not in use.
- Hand wash in warm water, towel off, and air dry in shade. Do not use any substances containing chlorine.

AT THE END OF ITS USE LIFE, DISCARD IT ALONG WITH YOUR NORMAL WASTE

INDICATIONS

- Degenerative, Herniated or Bulging Disc
- Spondylolysis
- Spondylolisthesis
- Facet Syndrome
- Spinal Stenosis
- Tumors or Inflammation
- Chronic back pain
- Spinal Instability
- Lumbago

And other conditions determined by an orthopedic surgeon, orthotist or other healthcare professional

Subscribe to our Fitting Instructions Videos on YouTube at:
www.youtube.com/OPTECUSA

OPTEC USA

975 Progress Circle Lawrenceville GA 30043
Ph: 1.770.513.7380 • F: 1.770.513.4365
Ph: 888.982.8181 • F: 877.288.4197
Email: sales@optecusa.com



www.optecusa.com

PRO & X



FITTING INSTRUCTIONS



VENUM 1
PRO LSO



VENUM 1
X LSO

