

Premier Series Locking

Specialty Bracing Solutions

Contact Information

☐ Clinician ☐ Fitter/Assistant/Tech ☐ Other: _____

Name: _____

Email: _____ Phone: _____

Ordering Clinician

☐ CPO ☐ CO ☐ CP ☐ Other: _____

Name: _____

Email: _____ Phone: _____

Billing & Shipping

PO#: _____

Billing Account#: _____

Shipping Account#: _____

Shipping Address: _____

City: _____ State: _____ Zip: _____

Shipping Preference

☐ Ground

☐ Next Day AM

☐ Next Day PM

☐ 2-Day AM

☐ 2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

Patient Information

Fit Date: _____ Patient ID: _____

Age _____ Weight _____ (LBS) Height _____ (IN)

☐ Male ☐ Female Leg: ☐ Left ☐ Right

Diagnosis: _____

Surgeries (type/date): _____

Please Follow Step-By-Step Cast Protocol Instructions

M-L measurement at knee center _____

Casting Protocol: 18-20 inch length; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

Model

☐ Anterior Band Model



☐ Polio Band Model



Thigh Band Length ☐ 7" ☐ 8" ☐ 9" ☐ Other _____"

Tibia Band Length ☐ 7" ☐ 8" ☐ 9" ☐ Other _____"

☐ Single Strut KAFO with Heel Cup
(Must complete graphite extension to floor form)

☐ C/S Package* Only available on option 1
(for dynamic compression and enhanced suspension)

Hinge Options

☐ Set hinges to LOCK at cast position

_____ OR _____

☐ Set hinge at:

☐ 0° ☐ 5° ☐ 10° ☐ 15° ☐ 5° ☐ 10° ☐ 15°

☐ Other _____°

Townsend Hinge U0013

☐ 5 Bar Locking Joint

☐ Single Pivot Locking

Lock options

☐ Pull Up Lever

☐ Lateral side (recommended)

☐ Centered on Top Band

☐ Manual Trigger

☐ Free motion (5 bar locking only)

☐ No free motion

☐ Cable with twist release

☐ Positioned on anterior of shell

☐ Lateral side

(must select if posterior thigh shell)

☐ Free motion

☐ No free motion

Becker Hinge U0013-B

☐ Modular Ring Lock Model 1402-B

☐ Automatic Angled Levered Lock Model 1017A

☐ Modular Ratchet Lock Model 1018A

☐ Bend Levers As A Bail Rod

Lock Release Option

☐ BLISS

☐ Bail Rod

Townsend Joints Accessory

☐ Install Extension assist Bands/Posts

☐ Extension Stop Kit (5 Bar Free & TM5 only)

☐ Flexion Stop Kit (TM5 only)

☐ Flexion Stop for 5 Bar Free and 5 Bar Locking (Factory installed only)

☐ 15° ☐ 30° ☐ 45° ☐ 60° ☐ 75° ☐ 90°

Brace Color (Select one)

Textured Powdercoat Finish (Lightest, Most Durable Finish)

☐ Black

☐ Royal Blue

☐ Antique Pewter (Silver)

☐ Burgundy

Paint Finish

☐ High Gloss

☐ Beige

☐ Dark Violet

☐ Black

☐ Emerald Green

☐ Steel Blue

☐ Royal Blue

☐ White

☐ Indy Yellow

☐ Burgundy

☐ Burnt Orange

☐ Quicksilver

Custom High Gloss Paint Finish*

☐ Provide Custom Paint # _____

Options

Condylar pads: ☐ No ☐ Medial ☐ Lateral ☐ Both

☐ Anti-Migration Silicon Infused Strap Pads*

Sleeves* ☐ 18" Cotton ☐ 18" Neoprene ☐ 22" Neoprene

☐ 1/6" Comfort Thigh Sleeve

☐ C/S Wrap

(for compression & enhanced suspension)

*Indicates additional charges apply

OF-018 REV D

Received Date Thuasne USA's shipping department use only

Distributed by Thuasne USA
4615 Shepard Street, Bakersfield, CA 93313
Tele: 800.432.3466 • Fax: 800.798.2722
ThuasneUSA.com